

Consistency in Excellence.

CPG™ Gel Breast Implant Family Cohesive III™



CPG™ Breast Implant 331



CPG™ Breast Implant 332



CPG™ Breast Implant 333



CPG™ Breast Implant 321



CPG™ Breast Implant 322



CPG™ Breast Implant 323



CPG™ Breast Implant 311



CPG™ Breast Implant 312



CPG™ Breast Implant 313

CPG™ Gel Breast Implants

Simplicity itself.

CPG™ BREAST IMPLANT NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate Plus 2
	Low 1	Moderate 1

KEY: DIGIT 1 = COHESIVE LEVEL
DIGIT 2 = HEIGHT OF IMPLANT
DIGIT 3 = PROJECTION OF IMPLANT

Introducing a simple way to identify exactly what cohesive level, height and projection each implant provides. Three numbers are all you need to know.

The first number identifies the implant type. CPG™ Gel Breast Implants are Cohesive III™ and therefore the first number 3 stands for the cohesion level. The second number relates to the height of the implant, with 1 being the lowest height and 3 being the tallest.

The third number identifies implant projection, with 1 representing Moderate, 2 representing Moderate Plus and 3 representing High Projection. Ordering is as simple as 1 - 2 - 3.

A smarter approach to implant selection.

Implant selection should be as simple as it is accurate. With that goal in mind, Mentor developed the CPG™ Gel Breast Implant Family. Three heights, three projections and multiple ways to match your patient's unique physiology.

Designed to fit.

CPG™ Implants combine a streamlined inventory with highly flexible sizing. Essentially, physicians require less stock to match more body types.

- Designed to help create the natural shape of the breast
- Unique, tapered appearance / shape, with a gradually sloping upper pole and fuller lower pole
- MemoryGel™ Silicone Gel, Cohesive III™, our most cohesive, form stabilizing gel
- Highly compliant shell and orientation marks for ease of placement
- Implant shell designed to minimize silicone gel bleed
- SILTEX™ Texture to provide a disruptive surface for collagen interface
- Supported by a lifetime replacement policy. Additionally, it is supported by PatientSafe Coverage limited warranty or PatientCare Promise limited warranty*

*In case of eligibility, depending on the date of surgery, the implant will fall under one of the above warranties.

MemoryGel™ Breast Implants. All Cohesive, All The time™.

Mentor's proprietary cohesive silicone gel holds shape uniformly and has the natural 'give' of breast tissue. Beautifully pliant. Firm without rigidity. MemoryGel™ Implants offer patients a more natural look and feel. By varying gel components, Mentor produces a wide selection of consistencies from soft (Cohesive I™) to firm (Cohesive III™).



A cross-section of a CPG™ Gel Breast Implant, Cohesive III™

Important Safety Information:

MENTOR® MemoryGel™ Breast Implants and MENTOR® Saline-filled Breast Implants are indicated for breast augmentation in women (at least 18 years old) or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body, with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions, or are pregnant or nursing.

Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery. The most common complications with the MemoryGel™ Implant includes reoperation, capsular contracture, asymmetry, and breast pain. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.

The most common complications with MENTOR® Saline-filled Implants include reoperation, implant removal, capsular contracture, wrinkling, breast pain and deflation.

Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.


MENTOR® CONTOUR PROFILE™ Tissue Expanders and Smooth Round Tissue Expanders with Remote Dome are used for breast reconstruction following mastectomy. These expanders are intended for temporary subcutaneous or submuscular implantation and are not intended for use beyond six months. Do not use the CONTOUR PROFILE™ Tissue Expander and Smooth Tissue Expander with Remote Dome in patients where an MRI may be needed. The device could be moved by the MRI causing pain or displacement, potentially resulting in a revision surgery. The incidence of extrusion of the expander has been shown to increase when the expander has been placed in injured areas.

For detailed indications, contraindications, warning and precautions associated with the use of all Mentor implantable devices, which include MENTOR® Saline-filled Breast Implants, MENTOR® MemoryGel™ Breast Implants and Tissue Expanders, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorwllc.com.

CPG™ Gel Breast Implant Family

Featuring Cohesive III™ MemoryGel™ Silicone Gel

CPG™ BREAST IMPLANT NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate Plus 2
	Low 1	Moderate 1
KEY: DIGIT 1 = COHESIVE LEVEL DIGIT 2 = HEIGHT OF IMPLANT DIGIT 3 = PROJECTION OF IMPLANT		



CPG™ Breast Implant 331

Cohesive III™, Tall Height, Moderate Projection

The Tall Height, Moderate projection CPG™ Implant is designed for tall patients with a narrow torso seeking for moderate results in their augmentation procedure. This implant also shows excellent results when correcting tuberous breasts in the above mentioned patient type.

As a reconstruction device it can be used as a follow up implant after the Tall Height CONTOUR PROFILE™ Tissue Expander.

CPG™ Breast Implant 321

Cohesive III™, Medium Height, Moderate Projection

The Medium Height, Moderate projection CPG™ Implant is an excellent choice for most augmentation patients whether the augmentation procedure is associated with mastopexy or not.

It also provides an option in reconstruction patients when a high projection is not desired and can be used as a follow up implant after our Medium Height CONTOUR PROFILE™ Tissue Expander

CPG™ Breast Implant 311

Cohesive III™, Low Height Moderate Projection

The Low Height, Moderate projection CPG™ Implant is designed for augmentation or reconstruction procedures in short patients with a wide torso seeking regular projection.

As a reconstruction device it can be used after both our Low Height and Medium Height CONTOUR PROFILE™ Tissue Expander.

CPG™ Breast Implant 332

Cohesive III™, Tall Height, Moderate Plus Projection

The Tall Height, Moderate Plus projection CPG™ Implant has the same footprint as the CPG™ Implant 331 but with more inferior pole projection. It is an ideal device for patients with taller, narrow torsos desiring more projection.

This implant is also selected for reconstruction patients after using the Tall Height CONTOUR PROFILE™ Tissue Expander

CPG™ Breast Implant 322

Cohesive III™, Medium Height, Moderate Plus Projection

The Medium Height, Moderate Plus projection CPG™ Implant is suitable for a great number of augmentation patients. It offers Moderate Plus projection especially in the inferior pole.

As a reconstruction device it can be used after our Medium Height CONTOUR PROFILE™ Tissue Expander or associated with a Lateral Dorsi flap.

CPG™ Breast Implant 312

Cohesive III™, Low Height Moderate Plus Projection

The Low Height, Moderate Plus projection CPG™ Implant with greater inferior pole projection offers multiple options in breast surgery. It is a highly suitable implant for patients with a shorter, wider torso seeking breast augmentation with Moderate Plus projection especially in the lower pole.

As a reconstruction device it can be used after both our Low Height and Medium Height CONTOUR PROFILE™ Tissue Expander.

CPG™ Breast Implant 333

Cohesive III™, Tall Height, High Projection

The Tall Height, High projection CPG™ Implant provides tall height with the most projection. It is designed for augmentation procedures in tall, slender patients with a narrow torso when a high projection outcome is desired.

In a two stage reconstruction procedure the implant can be used after the Tall Height CONTOUR PROFILE™ Tissue Expander It is the ideal choice as well when having to match a large contralateral breast in single stage reconstruction cases.

CPG™ Breast Implant 323

Cohesive III™, Medium Height, High Projection

The Medium Height, High projection CPG™ Implant offers multiple options in augmentation procedures especially when correcting breast ptosis and no mastopexy is associated.

As well, it is an ideal device for reconstruction procedures where a high inferior pole projection is important. In a two stage reconstruction procedure the implant can be used after our Medium Height CONTOUR PROFILE™ Tissue Expander.

CPG™ Breast Implant 313

Cohesive III™, Low Height High Projection

The Low Height, High projection CPG™ Implant is an ideal augmentation device for short patients with a wide torso seeking high inferior pole projection.

In a two stage reconstruction procedure the implant can be used after both our Low Height and Medium Height CONTOUR PROFILE™ Tissue Expander.

