

TLC® Retractors

Self-Retaining Retractor Systems

Call or Email Today
800 869 7382
CS@AppliedMedical.net

TLC® Retractors

Self-Retaining Retractor System

The TLC® Self-Retaining Retractor System allows surgeons to retract and secure delicate soft tissue away from the operative surface, enabling greater focus on the procedure. TLC® retraction is continuous and lower profile when compared to manual retraction.




Faster Deployment

Hinges adjust to the patient's anatomy in an instant with the malleable function of the TLC® Retractors. There are no thumbscrews to adjust or obstruct the procedure.



Eliminate Slippage

TLC® ribbed stays provide a secure fit, virtually eliminating slippage and risk of popping out.



Increase Control

Uniform tension allows surgeon to better control the tissue as the retractor is precisely positioned around the surgical site.



Lower Profile

Increases control and exposure of even the smallest surgical sites.

TLC® Make Connections



TLC® Self-Retaining Retractor Systems

Retractor Stays	
Order #	Description
TLC103	Sharp Tip (5mm) Universal
TLC108	Blunt Tip (12mm) Universal
TLC150	Blunt Tip (12mm)
TLC110	Sharp Tip (5mm)
TLC100	Semi-Sharp Tip (5mm)
TLC090	Sharp Tip (3mm)
*TLC203	(20mm) 3-Finger Rake

Retractor Frames	
Order #	Description
TLC200	Large Size
TLC300	Large Size
TLC400	Med-Size
TLC500	Small Size
TLC600	Med-Size
TLC700	Small Size
TLC800	Small Size
TLC900	Large Size
TLC920	**Pediatric Size


Stays available in 4, 8 or 50 packs
*Stays only available in 4 packs

See frames below

**Can be used for pediatric procedures
Universal Stay 5mm and 12mm are recommended


TLC® Retractors

Wide range of self-retaining retractor systems



Disposable retractor rings are pre-sterilized and packaged singly in Tyvek® pouches for easy delivery to the sterile surgical field.

TLC® Retractor Stays




Stays available in 4, 8 or *50 packs

Stays are pouched sterile
4 and 8 packs come in one pouch
*50 packs are individually pouched

Controlled Versatility

The TLC® Retractors and Stays come in a variety of shapes and sizes to cover a broad spectrum of surgical procedures. This versatility will give the surgical assistant freedom to do other tasks, broadening the range of procedures the surgeon can perform at one time.



Specialized Surgical Kits Also Available

Wilson™ Enhanced Retractor System Kit
TLC5011 (Left)

Small Deaver – TLC601
17mm wide (Above)

Innovating. Educating. Changing Lives.™

Please call or email today for more information:

800 869 7382
CS@AppliedMedical.net

