

Tomorrow's Innovation

TITAN X450



The TITAN's 400 watts of Bright White light is the only lightsource on the market that delivers this amount of light to the surgical site.

- New Ceramic body elliptical lamp
- New Touch Screen Controls
- 1,000 hour lamp
- Dual Fan Cooling System
- External Air Filter
- Automatic Shutter
- Standby Mode
- New Mechanical Lamp Alignment
- Consistent High Lux Output
- Internal Self Diagnosis

TITAN X350

The TITAN's 300 watts of Bright White light is superior to any other 300 watt system in the market due to the advanced optic design.

- New Ceramic body elliptical lamp
- New Touch Screen Controls
- 1,000 hour lamp
- Dual Fan Cooling System
- External Air Filter
- Automatic Shutter
- Standby Mode
- New Mechanical Lamp Alignment
- Consistent High Lux Output
- Internal Self Diagnosis



TITAN X300



300 watts of Bright White light
Enhanced Quality
Whisper Quiet Operation

- Consistent High Lux Output
- 1,000 hour lamp
- Manual Intensity Slide
- Advanced Air Flow
- Automatic Shutter
- Reliable
- Whisper Quiet
- Quality, Performance and Advanced Aesthetics

TITAN X450 & X350

- Lamp Type Ceramic Type Xenon Lamp
- Lamp Power Rating 400 Watts & 300 Watts
- Lamp Color Temperature 5,600K
- Lamp Warranty 1000 Hours
- Front Panel Controls On/Off Switch, Touch Screen Controls, Turret
- Input Voltage Range 100 – 120 VAC 50 /60Hz, 220–240VAC 50Hz
- Power Consumption 720 Watts & 450 Watts
- Regulatory UL60601, EN 60601-1, EN 60601-1-2, CAN/CSA, C22.2, No 601.1-M90, CE
- Equipment Class BF-TYPE
- Dimensions 338mm (13.3") W x 155mm (6.1") H x 457mm (18") D
- Weight 10.0 kg (22 lbs)

TITAN X300

- Lamp Type Ceramic Type Xenon Lamp
- Lamp Power Rating 300 Watts
- Lamp Color Temperature 5,600K
- Lamp Warranty 1000 Hours
- Front Panel Controls On/Off Switch, Attenuator Control, Turret
- Input Voltage Range 100 – 120 VAC 50 /60Hz, 220–240VAC 50Hz
- Power Consumption 450 Watts
- Regulatory UL60601, EN 60601-1, EN 60601-1-2, CAN/CSA, C22.2, No 601.1-M90, CE
- Equipment Class BF-TYPE
- Dimensions 338mm (13.3") W x 155mm (6.1") H x 457mm (18") D
- Weight 10.0 kg (22 lbs)

Ordering Information

- TITANX450ST-K 400 Watt Xenon Lightsource, 4 Port Turret
 - TITANX350ST-K 300 Watt Xenon Lightsource, 4 Port Turret
 - TITANX300ST-K 300 Watt Xenon Lightsource, 4 Port Turret
 - SSL0070T Floorstand for TITAN Lightsources
 - SSX0450 Replacement Lamp Module for TITANX450
 - SSX0350 Replacement Lamp Module for TITANX350
 - SSX0300 Replacement Lamp Module for TITANX300
- (For International use replace -K with required power cord plug)



SUNOPTIC
SURGICAL®

Forhandles av:

Medistim Norge AS. Postboks 6631 Etterstad, 0607 Oslo
Tlf: 23 03 52 50. www.medistim.no