



Smooth Round Tissue Expanders




Device Volume cc	Diameter cm	Projection cm	Catalog #
400	11,3	5,9	350-4305M
550	13,4	5,8	350-4307M
700	14,6	6,7	350-4304M
850	15,0	7,2	350-4308M
1000	15,5	7,8	350-4306M

Smooth Crescent Tissue Expanders




Device Volume cc	Length cm	Width cm	Projection cm	Catalog #
50	9,0	5,7	2,7	350-4312M
100	11,0	6,9	3,7	350-4313M
200	14,6	7,5	4,3	350-4314M
300	16,5	10,3	4,7	350-4315M
400	18,2	11,4	5,2	350-4316M
500	19,6	12,4	5,5	350-4317M

Smooth Elliptical Tissue Expanders



Device Volume cc	Length cm	Width cm	Nominal Projection cm	TEMPORARY OVEREXPANSION		Catalog #
				Proj. cm	Max. Volume cc	
25	6,0	4,0	2,3	4,5	75	350-5301M
50	8,0	5,0	2,5	5,3	150	350-5302M
75	10,0	6,0	2,1	6,0	225	350-5303M
150	12,0	7,0	3,2	7,2	450	350-5304M
250	14,0	8,0	4,5	9,0	750	350-5305M
125	15,0	5,0	3,3	6,5	375	350-5306M

Smooth Rectangle Tissue Expanders



Device Volume cc	Length cm	Width cm	Projection cm	Catalog #
50	6,9	3,2	3,2	350-4309M
100	6,6	5,1	5,6	350-4301M
250	10,4	7,3	5,6	350-4302M
400	10,6	9,3	6,6	350-4303M
550	13,7	7,1	7,7	350-4310M
700	15,0	8,0	7,9	350-4311M