



Excellence through research

mondomed®

HIGH TECH MEDICAL PVA FOAM



MONDOCEL® ENT

**Supporting medical
professionals in taking care
of people** since 1983





PIONEERING KNOW-HOW IN PVA FOAM

MONDOMED® has been specializing in the production of high quality **polyvinyl alcohol (PVA) foam** products for medical applications and has more than 40 years of experience.

MONDOMED®'s continuous research and attention to quality has led to excellence and reliability.

Today, **MONDOMED®** is considered to be the only PVA producer in the world offering such a wide range of high tech PVA foam products for medical applications.

MONDOMED® PVA Benefits

POLYVINYL ALCOHOL (PVA) has unique technical and physical properties which makes the material very suitable for a wide range of medical applications.

MONDOMED® PVA foam tampons:

- Have a spongy, porous structure with a small pore size
- Minimize tissue ingrowth
- Have an impressive fluid absorption and retaining capacity
- Quickly expand after hydration and adapt to the anatomical shape
- Have a smooth foam texture and are atraumatic
- Are dust, latex and fibre-free
- Are biocompatible



Mondocel[®] NASAL **MONDOMED**[®] PVA TAMPONS

INDICATIONS

- **NASAL TAMPONS**

- Epistaxis
- Postoperative

- **SINUS TAMPONS**

- Postoperative

BENEFITS

- Very smooth, spongy porous structure
- Highly absorbent
- Highly compressed
- Quick expansion when hydrated
- Tamponade on the neighbouring tissues
- Haemostatic attributes
- Adapts well to anatomical shape
- Biocompatible
- Lint and fibre-free
- Small pore size minimizes tissue ingrowth
- Atraumatic removal
- A longitudinally inserted tube facilitates breathing



Mondocel[®]


NASAL

MONDOMED[®] PVA TAMPONS

AVAILABLE NASAL TAMPONS

STANDARD, WITH TUBE, WITH STRING, WITH TUBE & STRING

Article number	Dimensions (mm)	Tube	String	Packaging		
				Pcs/pouch	Pouches/box	Boxes/transport box
NASAL TAMPONS						
1155500	45x15x20	-	-	1	20	12
1155505	80x15x20	-	-			
1155500F	45x15x20	-	X			
1155505F	80x15x20	-	X			
1155510	45x15x20	X	-	1	20	8
1155515	80x15x20	X	-			
1155510F	45x15x20	X	X			
1155515F	80x15x20	X	X			



AVAILABLE PREFORMED NASAL TAMPONS

STANDARD, WITH TUBE, WITH STRING, WITH TUBE & STRING

Article number	Dimensions (mm)	Tube	String	Packaging		
				Pcs/pouch	Pouches/box	Boxes/transport box
PREFORMED NASAL TAMPONS						
1155518	80x15x30	-	-	1	20	12
1155519	35x15x25	-	-			
1155518F	80x15x30	-	X			
1155519F	35x15x25	-	X			
1155516	80x15x30	X	-	1	20	8
1155516F	80x15x30	X	X			



Mondocel[®]

NASAL

MONDOMED[®] PVA TAMPONS

AVAILABLE SLIM PREFORMED NASAL TAMPONS

STANDARD, WITH STRING

Article number	Dimensions (mm)	String	Packaging		
			Pcs/pouch	Pouches/box	Boxes/transport box
SLIM PREFORMED NASAL TAMPONS					
1155522	80x10x30	-	1	20	12
1155522F	80x10x30	X			
SLIM PREFORMED NASAL TAMPONS (PEDIATRIC)					
1155517	55x10x25	-	1	20	12
1155517F	55x10x25	X			



AVAILABLE EPISTAXIS NASAL TAMPONS

STANDARD, WITH STRING

Article number	Dimensions (mm)	String	Packaging		
			Pcs/pouch	Pouches/box	Boxes/transport box
EPISTAXIS NASAL TAMPONS					
1155506	100x15x25	-	1	20	12
1155506F	100x15x25	X			
1155521	55x15x25	-			
1155521F	55x15x25	X			



Mondocel[®]
MONDOMED[®] PVA TAMPONS

SINUS

AVAILABLE SINUS TAMPONS

STANDARD, WITH STRING

Article number	Dimensions (mm)	String	Packaging		
			Pcs/pouch	Pouches/box	Boxes/transport box
SINUS TAMPONS					
1155524F	35x12x20	X	2	20	12
1155525F	25x12x20	X			
1155526F	35x9x20	X			
1155527F	25x9x20	X			
1155528F	45x9x12	X			
1155529F	45x6x12	X			
1155530F	35x9x12	X			
1155531F	35x6x12	X			





Mondocel[®] EAR MONDOMED[®] PVA TAMPONS

INDICATIONS

- **EAR TAMPONS**
 - Otitis externa
 - Postoperative after tympanoplasty

BENEFITS

- Very smooth, spongy porous structure
- Can be easily inserted into the auditory meatus
- Highly compressed
- Quick expansion when hydrated (highly absorbent)
- Tamponade on the neighbouring tissues
- Haemostatic attributes
- Adapts well to anatomical shape
- Biocompatible
- Lint and fibre-free
- Small pore size minimizes tissue ingrowth
- Atraumatic removal



AVAILABLE EAR TAMPONS

STANDARD, FENESTRATED, WITH STRING, FENESTRATED & STRING

Article number	Dimensions (mm)	Fenestrated	String	Packaging		
				Pcs/pouch	Pouches/box	Boxes/transport box
PREFORMED EAR TAMPONS						
1155555	15x12	-	-	1	20	12
1155556	13x9	-	-			
1155557	13x9	X	-			
1155579	15x9	X	-	1	50	8
1155580	15x9	-	-			
1155585	24x9	-	-	1	20	12
1155586	24x12	-	-			
1155555F	15x12	-	X			
1155556F	13x9	-	X			
1155557F	13x9	X	X			

MONDOCEL[®] ENT

**Supporting medical
professionals in taking care
of people** since 1983

The information, pictograms, photographs, graphics, texts, etc. contained in this catalogue are protected by industrial and/or intellectual rights and are the property of Mondomed nv.

As such, any reproduction, representation, adaptation, translation and/or modification in whole or in part, is prohibited.

© 2026 All rights reserved MONDOMED[®]

Produced in Belgium

Mondomed nv

Middenweg 20 3930, Hamont-Achel, Belgium

+32 11 44 09 90 | sales@mondomed.be | www.mondomed.be