

MONDOCEL®

Operating Room DEVICES



SURGICAL HEADBAND

● Made from high tech MONDOMED® PVA, the MONDOCEL® surgical headband is fiber free, hypoallergenic and highly absorbent. It will help the staff in the O.R. by quickly absorbing sweat and producing a cooling effect.

◆ Ref.: 1155450

Product 1155450 is packed as follows:

1 piece/package - 20 packages/box - 9 boxes/transportbox



DIAMOND KNIFE CLEANING BLOCK

● The MONDOCEL® PVA Diamond Knife Cleaning Block is fiber-free and has a high absorption capacity. It helps to remove blood and cellular debris to eliminate damage to the diamond knife during the cleaning process.

◆ Ref.: 1155301

Product 1155301 is packed as follows:

1 piece/package - 20 packages/box - 8 boxes/transportbox



MONDOMED® PVA foam devices are medical products with convincing technical properties:

- Fiber-free
- Highly absorbent
- Biocompatible
- High in purity
- High tensile strength
- Strong/non-shredding
- Environmentally friendly (biodegradable)




CLEAN PAD

● The MONDOMED® PVA Clean Pad is used to easily and efficiently remove fluids and debris from otolaryngological instruments. The MONDOCEL® Clean Pad is a highly absorbent instrument wiper, is fiber free and doesn't damage delicate instruments.

◆ Ref.: 1155280

Product 1155280 is packed as follows:

1 piece/package - 20 packages/box - 8 boxes/transportbox



MONDOMED nv
Middenweg 20
3930 HAMONT-ACHEL
BELGIUM

SUPPORTING MEDICAL PROFESSIONALS IN TAKING CARE OF PEOPLE... since 1983

WWW.MONDOMED.BE