

PRODUCT SPECIFICATION

Product name	Non-Woven Eye Dressing Pad
Product Ref	1152
Marking	Product name, Code Ref, CE mark, Lot number, Expiry Date, latex free, sterile and single use symbols.
Dimension	6x8cm
Material	Non-Woven / Viscose, heat sealed edges



Primary use	A single use, sterile eye dressing pad intended to protect the eye or absorb eye secretions.
--------------------	--

PACKAGING

MOQ	Inner Box Configuration
400	8 inner boxes of 50pcs

Quality certification	ISO 13485
Conformity to Regulations	Class Is – Medical Device Directive 93/42/EEC
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	5 Years if stored as recommended above
Additional information	This product is Latex Free and sterile.

Medical Device and Reporting System:

Incidents caused by this device must be reported immediately to our offices.

REVISION:	1.0	DATE	25/09/2020	AUTHOR	A Gidley
------------------	-----	-------------	------------	---------------	----------

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare T: +44 (0) 1530 278 111 W: www.365healthcare.com	Unit 6, Delta Park Industrial Estate Millmarsh Lane, London Enfield, EN3 7QJ
--	--

365 Healthcare is a trading division of Bunzl Retail & Healthcare Supplies
Limited Registered Office: York House, 45 Seymour Street, London W1H 7JT. Registered in England No. 6253