

## MEDILOOPS® SILICONE VESSEL LOOPS

### FIELDS OF APPLICATION:

- General Surgery
- Vessel Surgery
- Cardiac and Thoracic Surgery
- Neurosurgery
- Urology

### MATERIAL:

Vessel Loop: Silicone, radiopaque





### DIMENSIONS:

Length Vessel Loop: 50 cm





### PACKAGING:

double packed, sterile



### MEDILOOPS® SILICONE VESSEL LOOPS, MINI

| Ref.No. | Dimensions   | Application | Colour  | Packing Unit |
|---------|--------------|-------------|---|--------------|
| 01.2801 | 0.9 x 1.4 mm | veins       |    | 10 x 2 each  |
| 01.2802 | 0.9 x 1.4 mm | arteries    |    | 10 x 2 each  |
| 01.2803 | 0.9 x 1.4 mm | nerves      |  | 10 x 2 each  |
| 01.2804 | 0.9 x 1.4 mm | ureter      |  | 10 x 2 each  |

### MEDILOOPS® SILICONE VESSEL LOOPS, MAXI

| Ref.No. | Dimensions   | Application | Colour  | Packing Unit |
|---------|--------------|-------------|---|--------------|
| 01.2811 | 1.1 x 2.5 mm | veins       |  | 10 x 2 each  |
| 01.2812 | 1.1 x 2.5 mm | arteries    |  | 10 x 2 each  |
| 01.2813 | 1.1 x 2.5 mm | nerves      |  | 10 x 2 each  |
| 01.2814 | 1.1 x 2.5 mm | ureter      |  | 10 x 2 each  |

### MEDILOOPS® SILICONE VESSEL LOOPS, SUPERMAXI

| Ref.No. | Dimensions   | Application | Colour  | Packing Unit |
|---------|--------------|-------------|---|--------------|
| 01.2821 | 1.3 x 4.9 mm | veins       |  | 10 x 2 each  |
| 01.2822 | 1.3 x 4.9 mm | arteries    |  | 10 x 2 each  |