

# RingProtect<sup>TM</sup>

**Safe retraction, full protection**



# RingProtect™



RingProtect™ Disposable Retractor  
**Adjustable**



RingProtect™ Disposable Retractor  
**Non-adjustable**

## RingProtect™

### Key features

RingProtect™ Disposable Retractor is a great solution designed to retract surgical wound during surgery what exposes internal organs for operation and helps protect wound site from bacterial contamination or tumor cells implantation what provides substantial cost savings.

Circumferential elastic retraction maximizes the working area for better visualization and increased exposure.

RingProtect™ Disposable Retractor is perfect choice for all abdominal procedures as well as for thoracic and cardiac ones. Seamless and hassle-free conversion process when switching from competitive products to Grena RingProtect™ Disposable Retractor, as no clinical technique changes are required.

REF	Product	Rings Ø (mm)	Heigth (mm)	PCS per box
0221-060070150	Disposable Retractor 60x70x150 adjustable	60 x 70	150	5
0221-080090150	Disposable Retractor 80x90x150 adjustable	80 x 90	150	5
0221-120130150	Disposable Retractor 120x130x150 adjustable	120 x 130	150	5
0221-150160150	Disposable Retractor 150x160x150 adjustable	150 x 160	150	5
0221-180190200	Disposable Retractor 180x190x200 adjustable	180 x 190	200	5
0221-220230200	Disposable Retractor 220x230x200 adjustable	220 x 230	200	5
0221-270280250	Disposable Retractor 270x280x250 adjustable	270 x 280	250	5
0221-320330250	Disposable Retractor 320x330x250 adjustable	320 x 330	250	5
0221-040030025NA	Disposable Retractor 40x30x25 non-adjustable	40 x 30	25	5
0221-040030035NA	Disposable Retractor 40x30x35 non-adjustable	40 x 30	35	5
0221-050040025NA	Disposable Retractor 50x40x25 non-adjustable	50 x 40	25	5
0221-050040035NA	Disposable Retractor 50x40x35 non-adjustable	50 x 40	35	5
0221-070070035NA	Disposable Retractor 70x70x35 non-adjustable	70 x 70	35	5