

# UROMED »HYDROPUR®« STEUERBAR

## REF 4485 – 4488

### Zylindrische Spitze mit Perforation im Ureterteil

### Katheter zur inneren Schienung des Ureters


Für die ureterorenoskopische Anwendung ist die Ureter-Dauerschleife in Ch. 4,8 vorgesehen.

- aromatisches Polyurethan
- mit vergüteter Oberfläche
- geringe Inkrustationsneigung
- zur Langzeitanwendung geeignet
- gutes Vorschiebehaviorn durch mittlere Materialkonsistenz
- konische, dilatierende Spitze
- kleine, kreisrunde Pigtails
- atraumatische Augen
- vormontiertes Set zur einfachen Handhabung
- röntgenologisch darstellbar
- graduiert

### Catheter for the internal splinting of the ureter

The ureteral stent in Ch. 4,8 is destined for ureterorenoscopic application.

- aromatic polyurethane
- with refined surface
- low incrustation rate
- suitable for long-term application
- good advancement control because of medium material consistency
- conical, dilating tip
- small-sized, circular pigtails
- atraumatic eyes
- pre-assembled set for easy handling
- overall radiopaque
- graduated
- manoeuvrable



Ureter-Dauerschleife beidseitiger Pigtail  
*ureteral stent, double pigtail cylindrical tip, with perforation in the ureteral part*

Ch. 6-Ch. 8  
Pushersystem  
Länge 45 cm  
*Ch. 6-Ch. 8  
pusher system  
length 45 cm*

Ch. 4,8  
Pushersystem  
Länge 70 cm  
*Ch. 4,8  
pusher system  
length 70 cm*

Spiralmandrin  
PTFE-beschichtet  
*guidewire  
PTFE coated*

1 Fixierklemme bei REF 4485  
*1 fixing clamp for REF 4485*

REF <i>REF</i>	Dauerschleife <i>Stent</i>	Mandrin <i>Guidewire</i>	PZN <i>PZN</i>
4485 - 24	Ch. 4,8 / 24 cm	0,028" / 150 cm	11371739
4485 - 26	Ch. 4,8 / 26 cm	0,028" / 150 cm	11371745
4485 - 28	Ch. 4,8 / 28 cm	0,028" / 150 cm	11371751
4485 - 30	Ch. 4,8 / 30 cm	0,028" / 150 cm	11371768
4486 - 24	Ch. 6,0 / 24 cm	0,035" / 100 cm	11371774
4486 - 26	Ch. 6,0 / 26 cm	0,035" / 100 cm	11371780
4486 - 28	Ch. 6,0 / 28 cm	0,035" / 100 cm	11371797
4486 - 30	Ch. 6,0 / 30 cm	0,035" / 100 cm	11371805

REF <i>REF</i>	Dauerschleife <i>Stent</i>	Mandrin <i>Guidewire</i>	PZN <i>PZN</i>
4487 - 24	Ch. 7,0 / 24 cm	0,035" / 100 cm	11371811
4487 - 26	Ch. 7,0 / 26 cm	0,035" / 100 cm	11371828
4487 - 28	Ch. 7,0 / 28 cm	0,035" / 100 cm	11371834
4487 - 30	Ch. 7,0 / 30 cm	0,035" / 100 cm	11371840
4488 - 26	Ch. 8,0 / 26 cm	0,038" / 100 cm	11371857
4488 - 28	Ch. 8,0 / 28 cm	0,038" / 100 cm	11371863
4488 - 30	Ch. 8,0 / 30 cm	0,038" / 100 cm	11371886

**Verpackung:** doppelt, steril  
**Original-Karton:** 10 Stück

**Packing:** double, sterile  
**Original carton:** 10 pieces

**UROMED**  
PRODUKTE FÜR DIE UROLOGIE