



# TESTICULAR IMPLANTS

**SMOOTH SURFACE, NATURGEL 1**



Sebbin  
PARIS

The partner your body can trust



1

### INDICATIONS

Testicular implants are used for reconstruction, in cases of agenesis or testicular atrophy, or after testicular removal (orchidectomy). They are also indicated for gender reassignment surgery.<sup>(1)</sup>

2

### BENEFITS

**Egg-shaped**, to closely match the morphology of a testicle.

**Smooth-surfaced** envelope, pre-filled of Naturgel 1, a medical grade cohesive gel.

**A double occlusion patch with a removable thread** to guide the needle and prevent a puncture of the implant.

**A range featuring five sizes** (variable height, volume and diameter), corresponding to each patient's morphology.

Source:  
(1) Pigot GLS, Al-Tamimi M, Ronkes B, van der Sluis TM, Ozer M, Smit JM, et al. Surgical Outcomes of Neoscrotal Augmentation with Testicular Prostheses in Transgender Men. J Sex Med. 2019;16(10):1664-71.




## Testicular implants

Our range includes five different sizes of smooth-surfaced testicular implants to match the morphology of each patient:



Height (mm)	Diameter (mm)	Volume (ml)	Reference
25	21	6	LS 01 025
34	27	12	LS 01 034
40	30	17	LS 01 040
50	33	27	LS 01 050
53	38	38	LS 01 053

The dimensions of our implants may slightly vary from the data mentioned in this document.



The fixation guide is not implantable and must be removed after having fulfilled its function.



The partner your body can trust

**Sebbin Benelux SPRL**

Tél.: +32 2 792 57 76 - Fax: +32 2 792 57 73  
Email: info@sebbin.be

**Sebbin Deutschland GmbH**

Tél.: +49 2102 70 88 00-23 - Fax: +49 2102 70 88 00-99  
Email for Germany: service@sebbin.de  
Email for Austria: austria@sebbin.de

**Groupe Sebbin SAS, France**

Tél.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88  
Email: serviceclients@sebbin.com

**Sebbin Ibérica Productos Médicos, S.L.**

Tél.: +34 93 159 3003 - Fax: +34 93 676 3053  
Email: info@sebbin.es

**Sebbin Suisse Sarl**

Tél.: +41 43 336 72 30 - Fax: +41 43 336 72 36  
Email: info@sebbin.ch

**Sebbin UK Ltd**

Tél.: 0800 652 4028 - Fax: +4433 0088 9575  
Email: custserv@sebbin.uk

This document is intended for healthcare professionals. According to the Medical Device Directive 93/42/EEC, testicular implants are Class IIb medical devices. They are intended to be used in plastic, reconstructive and aesthetic surgery. Our testicular implants are sold by the unit and delivered sterile, for single-use. They are manufactured with medical grade raw materials, biocompatible and implantable in the long-term. Our manufacturing processes follow the NF EN ISO 9001: 2015 and NF EN ISO 13485: 2016 standards. Please read carefully the instructions for use of this device before using. They are available on [www.mysebbin.com](http://www.mysebbin.com).

[www.sebbin.com](http://www.sebbin.com)

GRUPE SEBBIN SAS

39 Parc d'Activités des Quatre Vents - CS80010 Boissy l'Aillierie - 95651 CERGY-PONTOISE CEDEX  
Tel.: + 33 1 34 42 13 28 - Fax: + 33 1 34 42 16 88