


LAP-5[®]

Disposable Endoscopic Instruments Scissors, Dissectors, Graspers

Precision and ergonomics in your hand



reddot winner 2024
medical devices and technology design



PROVIDING SPECIALISTS WITH
INNOVATIVE
HEALTHCARE SOLUTIONS

WITH THE CENTRAL FOCUS ON
PATIENTS' SAFETY AND USERS' COMFORT



LAP-5[®] Disposable Endoscopic Instruments

Broad spectrum of newly designed monopolar disposable endoscopic instruments was developed to enable minimally invasive surgery procedures in a variety of specialities.

The unique, modern and ergonomical single hand operation design facilitates performance of even the most complicated both standard as well as bariatric minimally invasive procedures.

Sharp scissors allow comfort cutting with coagulation function.



A wide range of forceps gives the surgeon pos-

sibility to choose the optimal instrument for a given procedure. Switchable ratchet contributes to a significant improvement in the comfort of the surgeon's work.

The narrow profile of the dissector jaws together with their wide opening angle enables effective separation of tissue structures.

Each instrument has a standard 4mm HF connector for coagulation during the procedure.

Product delivered sterile for single use.



Scissors



Dissectors









Graspers

LAP-5[®] is a completely redesigned concept of disposable endoscopic surgical instrument.



LAP-5[®] is assembled without any use of glue or chemicals to join the parts. Every component of the instrument is designed and assembled on a "CLICK" basis which removes the necessity to use any chemical-based adhesives. This also make LAP-5[®] less heavy on environment when it becomes hospital waste after the procedure.























Scissors

REF	PRODUCT	TIP	TYPE/ COLOUR	SIZE	PCS PER BOX
0208-DSLAP501	LAP-5® Disposable Endoscopic Metzenbaum Scissors		Endo ●	5 mm	4/16
0208-DSLAP502	LAP-5® Disposable Endoscopic Straight Scissors		Endo ●	5 mm	4/16
0208-DSLAP503	LAP-5® Disposable Endoscopic Hook Scissors		Endo ●	5 mm	4/16
0208-DSLAP501B	LAP-5® Disposable Endoscopic Metzenbaum Scissors		Endo 45 ●	5 mm	4/16
0208-DSLAP502B	LAP-5® Disposable Endoscopic Straight Scissors		Endo 45 ●	5 mm	4/16
0208-DSLAP503B	LAP-5® Disposable Endoscopic Hook Scissors		Endo 45 ●	5 mm	4/16

Dissectors

REF	PRODUCT	TIP	TYPE/ COLOUR	SIZE	PCS PER BOX
0208-DDLAP501	LAP-5® Disposable Endoscopic Dissector Maryland		Endo ●	5 mm	4/16
0208-DDLAP501B	LAP-5™ Disposable Endoscopic Dissector Maryland		Endo 45 ●	5 mm	4/16

Graspers

REF	PRODUCT	TIP	TYPE/ COLOUR	SIZE	PCS PER BOX
0208-DGLAP501	LAP-5® Disposable Endoscopic Atraumatic Fenestrated Grasper		Endo 	5 mm	4/16
0208-DGLAP502	LAP-5® Disposable Endoscopic Grasper Allis		Endo 	5 mm	4/16
0208-DGLAP503	LAP-5® Disposable Endoscopic Grasper Maxi Grip		Endo 	5 mm	4/16
0208-DGLAP504	LAP-5® Disposable Endoscopic Toothed Grasper		Endo 	5 mm	4/16
0208-DGLAP505	LAP-5® Disposable Endoscopic Grasper Babcock		Endo 	5 mm	4/16
0208-DGLAP501B	LAP-5® Disposable Endoscopic Atraumatic Fenestrated Grasper		Endo 45 	5 mm	4/16
0208-DGLAP502B	LAP-5® Disposable Endoscopic Grasper Allis		Endo 45 	5 mm	4/16
0208-DGLAP503B	LAP-5® Disposable Endoscopic Grasper Maxi Grip		Endo 45 	5 mm	4/16
0208-DGLAP504B	LAP-5® Disposable Endoscopic Toothed Grasper		Endo 45 	5 mm	4/16
0208-DGLAP505B	LAP-5® Disposable Endoscopic Grasper Babcock		Endo 45 	5 mm	4/16



Grena Ltd.

1000 Great West Road, Brentford,
London TW8 9HH, United Kingdom



KONMEX sp. z o.o.

Suite 6, 63A Bartycka, Street,
00-716 Warsaw, Poland

**CUSTOMER SERVICE CENTER
NOTTINGHAM OFFICE, UK**

☎ +44 7791 633 336

✉ sales@grena.co.uk

www.grena.co.uk

