



SUPERIOR MEDICAL CAMERAS

# proCam XS

When size matters

**novoCam**   
SHARE *your* VISION

## proCam XS

A miniature size camera with giant 4K features for loupe or light mounting. The proCam XS is half the size of the previous generation and boasts a Sony IMX 415 sensor for pristine video and still images, weighs only 14 grams, with manual focus to ensure 100% clarity when working loupe mounted.

- Medical 4K miniature loupe camera
- Only 14 grams
- Mounts to most loupes, headbands and treatment lights (optional attachments)
- Optional customized light adapter
- USB connection
- 4K / 1080p / 720p @30fps; 8.3 MP photos

## Benefits

- Easy mounting on loupes, headbands or treatment lights
- Instant hands-free operation by foot pedal → no treatment interruptions
- Quick documentation and sharing
- Simple USB installation and integration
- Modern case visualization and patient treatment plans
- Certified and secured products, EU MDR 2017/745 Class I
- Fast ROI by improvements in treatment acceptance rates
- Powerful tool for safe and repeatable education

**novoCam** 

**www.novocam.com**  
Novocam Medical Innovations Oy

novoCam cameras are manufactured in Finland according MDR 2017 /745, ISO13485 and IEC60601 requirements and registered as Class I Medical Device in the EU.